




## Improved $Mg^{2+}$ -activated RT-PCR with ES112 & ES113



Three different thermostable DNA polymerases were used to reverse transcribe an HIV transcript RNA template and subsequently amplify the cDNA in a coupled RT-PCR in the presence of either 3 mM  $Mg^{2+}$  or 3 mM  $Mn^{2+}$ . After 55 cycles of PCR, gel results demonstrate specific amplification products from RNA with Z05 in the presence of  $Mn^{2+}$ , but no specific product was observed when  $Mg^{2+}$  was used as the divalent metal ion activator. However, designer enzymes ES112 and ES113 produced specific amplification product with either  $Mg^{2+}$  or  $Mn^{2+}$  activation.




## Reduced RT Time Requirement for ES112 & ES113 in Mn<sup>2+</sup>

A 280 bp GAPDH RNA template was subjected to various RT incubation times and then amplified by PCR. In all cases PCR profiles were identical and the results were analyzed by kinetic PCR. The C<sub>T</sub> values of growth curves are plotted in the following chart:




Following a 30 min RT incubation time and Mg<sup>2+</sup> activation, the mutant enzymes ES112 and ES113 achieved RT activity similar to Mn<sup>2+</sup>-activated wild-type ZO5 DNA polymerase. With Mn<sup>2+</sup> activation, the mutant enzymes exhibited similar RT activity, but with much shorter RT incubation times (as low as 5 min). Even with no added RT incubation time there were only slight C<sub>T</sub> delays for Mn<sup>2+</sup>-activated mutant enzyme amplifications and initial PCR ramp times apparently are sufficient for the RT step to occur.

## Efficient RT-PCR at Decreased ES112 & ES113 Enzyme Concentrations



Enzyme concentration was titrated from 30 U down to 2.5 U per reaction for ZO5, ES112 and ES113. A significantly higher  $C_T$  value is observed with 2.5 U of ZO5 when compared to higher enzyme concentrations. The ES112 and ES113 perform relatively efficient RT-PCR with as little as 2.5 U of enzyme per 50  $\mu$ L reaction.

## Improved Low Copy Sensitivity with ES112 in $Mn^{2+}$ -activated RT-PCR



### ES112, $Mg^{2+}$ 10/32 Positives


Nominally 0.5 copies of HIV transcript RNA per reaction were amplified in 50  $\mu$ L RT-PCR amplifications optimized for  $Mg^{2+}$ -activated ES112,  $Mn^{2+}$ -activated ES112 or  $Mn^{2+}$ -

### ZOS, $Mn^{2+}$ 10/32 Positives

activated ZOS ("Gold Standard"). The  $T_m$  of end-point RT-PCR product was used to distinguish successful amplification of transcript RNA (specific product) from negative reactions (nonspecific product). The  $Mg^{2+}$ -activated ES112 reactions had the same low copy sensitivity as the  $Mn^{2+}$ -activated ZOS, while the low copy sensitivity was observed to be twice as good with  $Mn^{2+}$ -activated ES112.

### ES112, $Mn^{2+}$ 20/32 Positives

## RT-PCR Using $Mg^{2+}$ -activated CS6 DNA Polymerase



Various concentrations of pAV109 transcript RNA were amplified by single-buffer RT-PCR. All reactions contained 2 mM  $Mg^{2+}$  and CS6 DNA polymerase. Following 45 cycles of PCR, products of the correct size were observed with as little as  $10^3$  copies of RNA per reaction. Negative control reactions lacking RNA transcript produced no specific product of the expected amplicon size.